

# The Smart MR Injection System



Clear Direction.  From Diagnosis to Care.

**medrad**® MRXperion  
MR Injection System



# The Smart MR Injection System



Clear Direction.  From Diagnosis to Care.

**MEDRAD**® MRXperion  
MR Injection System

# Efficient. Smart.


The MEDRAD® MRXperion MR injection system has been designed to support you in providing efficient and accurate care.




- › Snap & Go syringe design
- › Auto Docking
- › Auto Loading
- › Auto Prime
- › Auto Retract

- › Convenience tray for syringe kit and supplies

- › Complete injection control within scan room
- › Illuminated fluid volume indicators



Sequential phase status



- › On-board eGFR calculators
- › Weight-based dosing calculator

# The Medrad MRXperion Workflow

## Standard Workflow

The new MEDRAD® MRXperion syringes with patented “snap & go” syringe design and the automated features simplify the preparation of the injector. It’s easy to use and saves time.

On board eGFR and weight-based dosing calculators provide important information that you may use to improve patient safety through individualized protocols via an intuitive, user friendly interface.

The test injection and KVO (keep vein open) function, together with injection pressure graphs, facilitate the high-quality performance of MEDRAD® MRXperion.

MR injection protocols can be stored for future reference ensuring consistent and repeatable results.



## Optimized Workflow with //Automated Documentation\*

**Smart ID (Barcode Reader)** captures contrast media data.

**Modality Worklist Import** reflects the full patient data for examination planning.

**PACS Interface** automatically generates a secondary capture, including injection and patient data.

**Dictation System Integration** automatically populates contrast and injection details.




\* Available as an option.

# Innovation in MR.

Optimize your workflow with //Automated Documentation\* and create new dimensions in digital healthcare.



Barcode reader\*\*



PACS Interface

## Benefit from improved quality of your data

- › Pull up patient data from the modality worklist
- › Provide the right injection protocol for each patient
- › Automatically capture contrast media information with the barcode reader

## Access data where you need it

- › Easily add all important injection and patient data to the patient's DICOM records
- › Simplify your reporting process, while eliminating potential manual errors

## Understand and analyze contrast usage

- › Perform statistical analysis
- › Centralize protocol management

## Maximize uptime with VirtualCare™\*

- › With the inclusion of a VirtualCare™ remote support system, alerts are monitored and resolved proactively

\* Available as an option.

\*\* Available in select markets.

# The Smart MR Injection System

Your Needs > Our Actions > Innovative Solutions

## Enabling Workflow-Efficiency<sup>1</sup>

From the way the syringes are installed to the way in which patient and injection data is made available to be utilized, the MEDRAD® MRXperion syringes and automated documentation\* solution have been designed to save you time.

**This is efficient.**

<sup>1</sup> using auto features of MRXperion



## Increases Confidence<sup>2</sup>

Reduce manual work and potential errors with automated documentation\*, while improving the accuracy and consistency of recorded protocols and associated data. MEDRAD® MRXperion gives you the confidence that you are acting on accurate information.

**This is smart.**

<sup>2</sup> dose is accurately delivered using MRXperion



## Helping You Focus On the Patient<sup>3</sup>

The combination of all the innovative features of MEDRAD® MRXperion means one thing – allowing you to provide the best care possible for your patient.

**This is MEDRAD® MRXperion.**



<sup>3</sup> auto features of MRXperion help minimize injector preparation time



\* Available as an option.

# Specifications

Components	Dimensions
<b>Control room unit</b>	15.12" (38.4 cm) W
	10.76" (27.33 cm) H
	7.94" (20.17 cm) D
<b>Scan room unit</b>	23.30" (59.0 cm) W
	71.40" (181.0 cm) H
	23.30" (59.0 cm) D
<b>Power supply</b>	7.60" (19.0 cm) W
	3.40" (9.0 cm) H
	15.40" (39.0 cm) D
<b>Electrical</b>	
<b>Voltage requirements</b>	100–240 VAC
	(50/60 Hz)
	120–210 VA
<b>Electrical leakage</b>	Unit <100 microamperes
	Patient <100 microamperes
	Earth <500 microamperes
<b>System capabilities</b>	
<b>Syringe capacities</b>	<b>Syringe A:</b> 65 mL
	<b>Syringe B:</b> 115 mL
<b>Programmable volume range (mL)</b>	<b>Syringe A:</b> 0.5 mL to max syringe volume in: 0.1 mL increments from 0.5 mL to 31 mL 1 mL increments above 31 mL
	<b>Syringe B:</b> 1 mL to max syringe volume in 1 mL increments
<b>Programmable flow range (mL/sec)</b>	0.01 to 10 mL/s in: 0.01 mL/s increments between 0.01 and 3.1 mL/s 0.1 mL/s increments between 3.1 and 10 mL/s
<b>KVO (keep vein open)</b>	6 factory presets of 0.25 mL every 15, 20, 30, 45, 60 or 75 sec
<b>Test inject</b>	Configurable from 0.5 mL to 20 mL in 0.1 mL increments
<b>Pressure range (psi)</b>	6 factory presets from 100 to 325 PSI or, 690 to 2240 kPa
<b>Injection/post injection reminders</b>	Up to 5 settings of 1 sec to 20 minutes in 1 sec increments
<b>Injection protocol storage</b>	60 protocols up to 6 phases each

<b>System capabilities (continued)</b>	
<b>Injection hold/pause</b>	Up to 20 minutes in 1 sec increments
<b>eGFR calculator</b>	<b>For adults:</b> MDRD, Cockcroft-Gault, Modified Cockcroft-Gault and CKD-EPI methods
	<b>For children:</b> Bedside Schwartz method
<b>Informatics connectivity capability</b>	Yes, with optional Radimetrics™ Enterprise Platform*
<b>Remote service capability</b>	Yes, with optional VirtualCare™ Remote Support*
<b>Mounting</b>	Floor pedestal – includes integrated IV pole
<b>Accessories</b>	
<b>IV Pole</b>	Choice of single hook or double hook IV pole
<b>Prefilled syringe adaptor</b>	Choice of 4 different adaptors
<b>Mobile clamp kit</b>	Available for mobile MRI scanners
<b>Ordering information</b>	
<b>Injector system</b>	Catalog # MRXP 200: MEDRAD® MRXperion MR Injection System
<b>Disposable syringe kit</b>	Catalog # XP 65/115 VS:
	<b>Dual syringe kit</b>
	<ul style="list-style-type: none"> <li>&gt; 65 mL Qwik-Fit syringe for contrast media</li> <li>&gt; 115 mL Qwik-Fit syringe for saline</li> <li>&gt; One large spike for saline</li> <li>&gt; One small spike for contrast media</li> <li>&gt; 96" low pressure connector tube with T-connector and check valve</li> </ul>
<b>Single syringe kit</b>	Catalog # XP 115 VS
	<b>Contents:</b> <ul style="list-style-type: none"> <li>&gt; 115 mL Qwik-Fit syringe for saline</li> <li>&gt; One large spike for saline</li> <li>&gt; 96" low pressure connector tube with T-connect or and check valve</li> </ul>

\* Radimetrics™ Enterprise Platform and VirtualCare™ Remote Support are options that can be ordered. They are not included as part of the base product.

**Clear Direction.  From Diagnosis to Care.**

Bayer reserves the right to modify the specifications and features described herein or to discontinue any product or service identified in this publication at any time without prior notice or obligation. Please contact your authorized Bayer representative for the most current information.

All patient data that appear in this document are fictitious. No actual patient information is shown.

Bayer, the Bayer Cross, MEDRAD MRXperion, MEDRAD, MRXperion, Radimetrics and VirtualCare are trademarks owned by and/or registered to Bayer in the U.S. and/or other countries. Other trademarks and company names mentioned herein are properties of their respective owners and are used herein solely for informational purposes. No relationship or endorsement should be inferred or implied.

© 2022 Bayer. This material may not be reproduced, displayed, modified, or distributed without the express prior written consent of Bayer.



More information on  
radiology.bayer.com



**Manufacturer**  
Bayer Medical Care Inc.  
1 Bayer Drive  
Indianola, PA 15051-0780  
U.S.A.  
Phone: +1-412-767-2400  
+1-800-633-7231  
Fax: +1-412-676-4120

Bayer Co. (Malaysia) Sdn. Bhd.  
Radiology, Pharmaceuticals  
25-03 & 25-04, Level 25, IMAZIUM  
No 8, Jalan SS21/37, Damansara Uptown  
47400 Petaling Jaya, Selangor, Malaysia.  
Phone: +603 7801 3088  
Fax: +603 7733 9272

**CE 2797**